

# NEXA AFFORDABILITY: MULTIPLE PATIENTS. ONE SOLUTION.



DFU

PI

VLU

## Smart Choices Add Up

When you're focused on effective, efficient patient care, the right partner can make a huge difference. NEXA is the affordable NPWT option for your facility.

### Here's why:



\*For device only. Additional costs for dressings and EcoPouch™ containers will vary depending on the number required for each patient.

## ORDERING Information

Product	Description	Part #	UOM
NEXA NPWT System Device Kit	NEXA therapy unit, power supply, carry case	414-00-300-GBR	1
NEXA Dressing -125mmHg	1 occlusive drape, 1 tubing connection pad, 1 foam filler (15cm x 7.5cm)	414-00-501	Box of 5
NEXA EcoPouch™ Fluid Container Pack	120ml flexible fluid container pouch and tubing	414-00-502	Box of 5
NEXA Drape	Occlusive drape individual sheets	414-00-600	Box of 10
NEXA Foam	Extra individual foam filler (15cm x 7.5cm)	414-00-700	Box of 5
Carrying Case	Carrying case for NEXA therapy unit	414-00-400	1
Power Supply	Power supply for NEXA therapy unit	414-00-451	1

## SIMPLE PORTABLE AFFORDABLE



Your patients deserve it all. Get the **NEXADVANTAGE**.

**Experience the Most Positive Thing to Happen to Negative Pressure Wound Therapy**



Visit [www.AOTInc.net](http://www.AOTInc.net) for more information.

AOTI is recognized as a leader in wound care technology around the world.



AOTI Ltd.  
Unit 20, Glenrock Business Park,  
Ballybane, Galway, Ireland.  
E: [ukcustomerservice@aotinc.net](mailto:ukcustomerservice@aotinc.net)  
P: +44 (808) 196-4364



# What is the **NEXA** ADVANTAGE?

The **NEXA Negative Pressure Wound Therapy (NPWT) System** has been developed to provide uncompromising negative pressure therapy for patients with chronic and acute wounds.

## ◆ **SIMPLE:**

- With no complicated interfaces or touchscreens, NEXA's ultra-simple, one-button operation requires minimal training and reduces the potential for user error
- Sleek, streamlined design with no maintenance required
- Accurate pressure control is built in, as the pre-set SmartFlow™ tubing connector constantly regulates negative pressure

## ◆ **PORTABLE:**

- Lightweight and small, ideal for mobile patients
- Easy to change, compact, lightweight pouch that collects fluid away from the wound
- Handy carrying case for patients on the go
- Battery lasts 10 hours before recharging

## ◆ **AFFORDABLE:**

- Designed to operate for up to 8 weeks (56 days or 1344 hours)
- Usable hours are depleted **only** when the device is providing therapy
- Remaining hours can be **used at your convenience on the next patient**
- **Immediate availability** with this off-the-shelf device that requires no tracking and incurs no runaway rental costs
- Disposable device eliminates servicing and maintenance costs

# A Powerful Balance of **SIMPLICITY** and **PERFORMANCE**

This **multi-week, multi-patient, disposable system** has three innovations that make all the difference.



### ◆ **Hydrodynamic Fluid Pump**

**Continues to remove fluid even in the presence of a minor dressing air leak**

- Constantly, actively removes exudate from wound, reducing potential for skin maceration
- No risk of internal device contamination because the disposable pump never touches the inside of the tubing or fluid container pouch



### ◆ **EcoPouch™**

**Draws exudate away from the wound site into the separate, disposable, flexible pouch**

- 120 ml capacity
- Includes fluid immobiliser
- Requires 58% less plastic content than conventional rigid canisters†



### ◆ **SmartFlow™ Connector**

**Pressure is constantly maintained at the wound**

- Maintains negative pressure levels preset to -125mmHg
- No need for complex pressure monitoring sensors or feedback loops within device, dressing, or tubing

