

**SIMPLE
PORTABLE
AFFORDABLE**



Your patients deserve it all. Get the **NEXADVANTAGE**.

**Experience the Most Positive Thing to Happen
to Negative Pressure Wound Therapy**



What is the **NEXADVANTAGE?**

The **NEXA™ Negative Pressure Wound Therapy (NPWT) System** has been developed to provide uncompromising negative pressure therapy for patients with chronic and acute wounds.

◆ **SIMPLE:**

- **With no complicated interfaces or touchscreens**, the ultra-simple, one-button operation requires minimal training and reduces the potential for user error
- **Sleek, streamlined design** with no maintenance required
- **Accurate pressure control is built in**, as the pre-set SmartFlow™ tubing connector constantly regulates negative pressure

◆ **PORTABLE:**

- **Lightweight and small**, ideal for mobile patients
- **Easy to change, compact, lightweight pouch** that collects fluid away from the wound
- **Handy carrying case** for patients on the go
- **Battery lasts 10 hours** before recharging

◆ **AFFORDABLE:**

- Designed to **operate for up to 8 weeks (56 days or 1344 hours)**
- Usable hours are depleted **only when the device is providing therapy**
- Remaining hours can be **used at your convenience on the next patient**
- **Immediate availability** with this off-the-shelf device that requires no tracking and incurs no runaway rental costs
- **Disposable device eliminates servicing** and maintenance costs

NEXAAFFORDABILITY: MULTIPLE PATIENTS. ONE SOLUTION.



Smart Choices Add Up

When you're focused on effective, efficient patient care, the right partner can make a huge difference. NEXA is the affordable NPWT option for your facility.

Here's why:



A Powerful Balance of **SIMPLICITY** and **PERFORMANCE**

This **multi-week, multi-patient, disposable system** has three innovations that make all the difference.



Hydrodynamic Fluid Pump

Continues to remove fluid even in the presence of a minor dressing air leak

- Constantly, actively removes exudate from wounds, reducing potential for skin maceration
- No risk of internal device contamination because the disposable pump never touches the inside of the tubing or fluid container pouch



EcoPouch™

Draws exudate away from the wound site into the separate, disposable, flexible pouch

- 120 ml capacity
- Much less shelf space required
- Includes fluid immobilizer
- Requires 58% less plastic content than conventional rigid canisters†



SmartFlow™ Connector

Pressure is constantly maintained at the wound

- Maintains negative pressure levels preset to -125mmHg
- No need for complex pressure monitoring sensors or feedback loops within device, dressing, or tubing



ORDERING Information

Product	Description	Part #	UOM
NEXA NPWT System Device Kit	NEXA therapy unit, power supply, carry case	414-00-300	1
NEXA Dressing -125mmHg	1 occlusive drape, 1 tubing connection pad, 1 foam filler (15cm x 7.5cm)	414-00-501	Box of 5
NEXA EcoPouch™ Fluid Container Pack	120ml flexible fluid container pouch and tubing	414-00-502	Box of 5
NEXA Drape	Occlusive drape individual sheets	414-00-600	Box of 10
Carrying Case	Carrying case for NEXA therapy unit	414-00-400	1
Power Supply	Power supply for NEXA therapy unit	414-00-450	1
Power Supply Extension	1.5 meter extension cable for power supply	414-02-023	1
NEXA Foam	Extra individual foam filler (15cm x 7.5cm)	414-00-700	Box of 5



3512 Seagate Way, Suite 100
Oceanside, CA 92056, USA
P: 1-760-431-4700 • F: 1-760-683-3063



Visit www.AOTInc.net
for more information.

AOTI is recognized as a leader in wound care
technology around the world.

